



**VASCULAR TECHNOLOGY
PROFESSIONAL PERFORMANCE GUIDELINES**

**Intracranial Cerebrovascular
Evaluations: Transcranial Doppler
Ultrasound and Transcranial Color
Duplex Imaging**

This Guideline was prepared by the Professional Guidelines Subcommittee of the Society for Vascular Ultrasound (SVU) as a template to aid the vascular technologist/sonographer and other interested parties. It implies a consensus of those substantially concerned with its scope and provisions. The guidelines contain recommendations only and should not be used as a sole basis to make medical practice decisions. This SVU Guideline may be revised or withdrawn at any time. The procedures of SVU require that action be taken to reaffirm, revise, or withdraw this Guideline no later than three years from the date of publication. Suggestions for improvement of this guideline are welcome and should be sent to the Executive Director of the Society for Vascular Ultrasound. No part of this guideline may be reproduced in any form, in an electronic retrieval system or otherwise, without the prior written permission of the publisher.

Sponsored and published by:
Society for Vascular Ultrasound
1601 Utica Ave S., Suite 213 Minneapolis, MN
55416
Tel.: 651-288-3431
Fax: 651-290-2266
E-mail: svuinfo@svunet.org
Internet: www.svunet.org

Copyright © by the Society for Vascular Ultrasound, 2026.
ALL RIGHTS RESERVED. PRINTED IN THE UNITED STATES OF AMERICA.
VASCULAR PROFESSIONAL PERFORMANCE GUIDELINE

Updated April 2026

PURPOSE

Transcranial Doppler ultrasound (TCD) and transcranial color duplex imaging (TCDI) studies measure blood flow velocities and direction within portions of the intracranial arteries. TCD and TCDI evaluate the anterior and posterior circulation territories to assess and manage patients with cerebrovascular disease. Testing is often performed to detect and monitor vasospasm after subarachnoid hemorrhage and screen children with sickle cell disease for stroke risk. Additionally, these exams can identify intracranial stenosis and occlusion. TCD can detect intracranial emboli and assess vasomotor reactivity.

APPROPRIATE INDICATIONS

Common indications for performance of this examination include, but are not limited to:

- Screening children with sickle cell disease for stroke risk
- Identification and monitoring of vasospasm after subarachnoid hemorrhage
- Detection of intracranial arterial stenosis or occlusion
- Assessment of cerebral vasomotor reactivity
- Support the clinical diagnosis of brain death
- Monitoring the effects of thrombolytic therapy in acute stroke patients
- Evaluation of extracranial stenosis on intracranial hemodynamics
- Detection of cerebral embolization
- Contrast enhanced detection of right-to-left cardiac/extracardiac shunts
- Intra-operative monitoring during carotid endarterectomy or coronary artery bypass graft procedures

CONTRAINDICATIONS AND LIMITATIONS

Contraindications and limitations may include the following:

- Inadequate acoustic windows and/or poor patient positioning
- Restless or uncooperative patients
- U.S. Food and Drug Administration (FDA) recommends the power output of the transducer be reduced to the lowest power that allows a complete orbital examination
- Remove contact lenses prior to the exam
- Output power must not exceed 10% of maximum emitted power or 17mW per cm² or equivalent
- In emergent or critical care situations, it may be difficult to perform a complete exam
- TCD/TCDI performance during interventional procedures requires additional precautions for the examiner, due to the presence of radiation
- Cerebral circulatory arrest evaluations should only be performed when skull is intact and the patient is normothermic, normocapnic and with BP above 100 mmHg systole. Patients with external ventricular drains will require clamping of the drain under supervision of the care team.

PATIENT COMMUNICATION

Prior to beginning the exam, the sonographer or examiner should:

- Introduce self and explain why the examination is being performed and indicate how much time the examination will take.
- Verify the patient's name and date of birth or utilize facility specific patient identifiers.
- Explain the procedure, taking into consideration the age and mental status of the patient and ensuring that the necessity for each portion of the evaluation is clearly understood.
- Respond to questions and concerns about any aspect of the evaluation.
- Educate the patient about risk factors for cerebrovascular disease (when appropriate).
- Refer specific diagnostic, treatment or prognosis questions to the patient's physician.
- Explain the necessity of and remove patient's eyeglasses, contacts and/or head-coverings.

PATIENT ASSESSMENT

Patient assessment should be completed before the evaluation is performed. This includes assessment of the patient's ability to tolerate the procedure and evaluation of any contraindications to the procedure. The sonographer or examiner should obtain a complete, pertinent history by interviewing the patient, the patient's representative, or a member of the care team overseeing management of the patient and review of the patient's medical record, if available. A pertinent history includes:

- Previous cardiovascular surgeries
- Current medications or therapies
- Presence of risk factors for cerebrovascular disease
- Presence of symptoms of cerebrovascular disease
- Results of other relevant diagnostic procedures
- Review of prior examinations to ensure the evaluation duplicates prior imaging and Doppler parameters
- Verify the requested procedure correlates with the patient's clinical presentation

PATIENT POSITIONING

Review patient positioning requirements for the examination and determine patient's ability to maintain proper positioning for all portions of the exam:

- Sitting or supine position with head supported or stabilized for transtemporal, orbital and sub-mandibular approach
- Turned to side with neck flexed/chin toward chest to optimize access to foramen magnum for sub-occipital (transforaminal) approach
- Alternatively, the sub-occipital exam can be performed with patient in sitting position, arms crossed and supported by stretcher or bedside table, head resting on arms so that neck is supported and relaxed
- Depending on neurological status, patients should be kept awake throughout

the exam due to changes in CO₂ during sleep which may affect velocities and result in misdiagnosis.

- The TCD examiner should be positioned at the patient's head with arms supported when possible.
- Inpatients may have restrictions on head positioning. A team member caring for an inpatient should be consulted prior to a TCD exam. Assistance may be required to maintain adequate head positions.

INSTRUMENTATION

Equipment used for TCD and TCDI is separate and distinct. The majority of velocity reference standards for intracranial studies were obtained at a zero-degree angle with non-imaging TCD equipment. Decreased velocities and inferior penetration are reported when comparing TCDI to TCD studies.

- ***Transcranial Doppler Ultrasound Equipment (TCD)***

- Utilizes a single crystal 1.5-2.5MHz pulsed Doppler probe
- Equipment and software specifically designed for TCD applications
- Displays all data in real time including:
 - Sample volume, sample depth
 - Time averaged maximum mean velocity
 - Peak systolic velocity/end diastolic velocity
 - Acceleration time
 - Pulsatility (PI) or resistive indices (RI)
 - Doppler power output and frequency
 - Direction sensitive Doppler blood flow meter
 - Bi-directional Doppler waveform display
 - Audible output and permanent recording of the spectral waveform
 - High intensity transient signals

- ***Transcranial Color Duplex Imaging Equipment (TCDI)***

- Duplex imaging scanner with B-Mode, color and spectral Doppler capabilities
- Utilizes a 1-5MHz phased-array transducer
 - Displays all data in real time including:
 - Sample volume, sample depth, Doppler angle
 - Time averaged maximum mean velocity
 - Peak systolic velocity /end diastolic velocity
 - Acceleration time

- Pulsatility (PI) or resistive indices (RI)
- Doppler power output and frequency
- Direction sensitive Doppler blood flow meter
- Bi-directional Doppler waveform display
- Visual display, audible output and permanent recording of the spectral waveforms with corresponding B-Mode/Color image

EXAM PROTOCOL

TCD studies follow a standard imaging protocol for the examination. A complete evaluation includes spectral Doppler analysis of all accessible portions of the major intracranial arteries and the distal cervical common carotid arteries. Velocities can be measured by automatic tracing or manual cursor placement. Bilateral evaluations are essential for a complete evaluation.

Transcranial Doppler- Complete Examination Protocol

TCD includes **Spectral Doppler waveforms** from the following vessels and *approaches* to obtain mean flow velocities (without angle correction):

Submandibular Window:

- Distal cervical internal carotid arteries (ICA)

Transtemporal Window:

- Terminal internal carotid arteries (TICA)
- M1 segment of the middle cerebral arteries (MCA)
 - Should obtain proximal, mid and distal measurements
- A1 segment of the anterior cerebral arteries (ACA)
- Anterior communicating artery, if detectable (ACoA)
- P1 and P2 segments of the posterior cerebral arteries
- Posterior communicating arteries, if detectable (PCoA)

Transforamenal/Suboccipital Window:

- Terminal vertebral arteries (VA)
- Proximal and distal segments of the basilar artery (BA)

Orbital Window: (may not be appropriate for all patients)

- Ophthalmic Artery (OA)
- Carotid siphon

Transcranial Doppler-Vasomotor Reactivity Protocol

Vasomotor reactivity testing is commonly ordered to evaluate cerebral reserve capacity in patients with unilateral carotid compromises and or basilar artery

compromises.

- A 1.5-2.5MHz transducer is used to insonate the MCA's either bilaterally or unilaterally. The patient is then asked to hold his breath for thirty seconds while the targeted artery is monitored.
- A normal response is for the mean flow velocity to reduce and then overshoot the baseline velocity as time expires and breathing resumes.
- Many stroke and TIA patients cannot hold their breath that long so to test reserve capacity an injection of Acetylzolimide (Diamox) can be used which takes about 10 minutes to maximally dilate the arteries, or a tank of air with high pCO₂ can be inhaled to precipitate vasodilation.
- This is the ultrasound version of a perfusion study without requiring a contrast medium.

Transcranial Doppler – Emboli Monitoring Protocol

Transcranial Doppler is the only modality that can detect and localize native emboli in real-time.

- A 1.5-2.5 blind pulsed-Doppler transducer is focused on the MCA and best done when fixed to a headband rather than handheld. Some headbands have servomotors inside them that allow for remote steering of the probe at the machine site. In 50-60% of stroke patients (reference?) bilateral MCA monitoring can be achieved. This can cut the exam time in half.
- Monitoring the basilar artery may be challenging in some individuals and is done with probe-in-hand. Because of this, ergonomic factors need to be considered when conducting a prolonged evaluation. Recording 500-750 beats per vessel may be adequate. Holding a pinch grip for more than 15 minutes is arduous.
- Detected artery to artery emboli will originate from an ipsilateral nidus.
- Emboli originating from a cardiac source may transit all the vessels targeted. M-mode spectral analysis can improve the display of the moving emboli as well as allow a calculation of emboli transit time.

Transcranial Doppler for Detection of Right-to-Left Cardiac Shunts

- One or both MCA's can be monitored with a 1.5-2.5 MHz pulsed-Doppler transducer.
- The transtemporal window should be evaluated for adequacy prior to placing an IV.
- Agitated saline should not be administered via PICC line unless specifically supported by departmental policy.
- Agitated saline administered via the jugular vein should be done by a qualified member of the patient's bedside care team.

- In the outpatient lab an IV line is placed by an institutionally qualified individual.
- On inpatient studies, the IV-line is often already placed but may be malfunctioning. Check for patency carefully and ask for assistance from the bedside care team if an IV appears to be malfunctioning.
- Different institutions have different guidelines for the ratio of air to saline in an agitated bolus. These may range from 9 cc saline with 1 cc of air to 9.8 cc of saline with 0.2 cc of air. Department protocol may restrict the number of boluses given to a patient and these may be determined by the total amount of air in delivered in each bolus. Pediatric patients may have further departmental restrictions. Recordings are made during the following conditions: at rest without Valsalva, at rest with Valsalva, after bolus injection, after bolus injection with Valsalva. The number of emboli are counted flowing through the target artery within 1 minute of the injection.
- A second injection is performed with a maximum Valsalva maneuver performed when the bolus is halfway injected. The number of microbubbles are counted flowing through the target arteries on this second maneuver.
 - Since the patient is awake and participating many shunts will be detected that are missed by transesophageal echocardiography.
- If the patient has no acoustic access through the temporal window monitoring the basilar artery is a suitable alternative.
- The transtemporal window should be evaluated for adequacy prior to placing an IV.

DISTINCTIONS BETWEEN TRANSCRANIAL DUPLEX IMAGING (TCDI) AND TCD EXAMS

TCDI combines real time B-Mode imaging of cerebral parenchyma with color Doppler imaging and spectral Doppler analysis of the major intracranial arteries. The standard examination protocol includes bilateral assessment of the distal extracranial ICA, ACA, MCA, PCA, vertebral and basilar arteries.

- Primarily used to assess stroke risk in children with sickle cell disease and to detect and monitor vasospasm following spontaneous subarachnoid hemorrhage
- Mean flow velocities and ratios remain the primary criteria for diagnosis of cerebrovascular disease
 - Waveforms obtained at 2-5mm increments with a 3-6mm sample volume size
 - Spectral analysis includes peak systolic velocity, end diastolic velocity, systolic upstroke or acceleration time, pulsatility index, and time-averaged maximum mean velocity
 - Angle corrected flow velocities may be higher than non-angle corrected flow velocities (should use a zero-degree angle)
 - Most diagnostic standards are based on non-imaging TCD studies

- Provides accurate B-Mode and color Doppler imaging of intracranial anatomy
 - Gray scale imaging of bony landmarks or parenchymal abnormalities
 - lesser wing of the sphenoid bone, petrous ridge of the temporal bone and cerebral peduncles
 - Color Doppler imaging of the major intracranial arteries
 - Should depict the ACA, MCA, PCA and vertebrobasilar arteries
 - May detect intracranial aneurysms $\geq 6\text{mm}$ in size
- Limitations compared to TCD remain:
 - poor beam penetration with lower vessel detection rates
 - suboptimal monitoring of emboli or vasomotor reactivity testing

Transcranial Doppler (TCD)

TCD provides spectral Doppler analysis of the major intracranial arteries. Diagnostic criteria should be based on published studies or internally validated criteria. The majority of reference standards were obtained with non-imaging TCD equipment. The standard examination protocol includes bilateral assessment of the distal extracranial ICA, ACA, MCA, PCA, vertebral and basilar arteries.

- TCD can be used to assess stroke risk in children with sickle cell disease and to detect and monitor vasospasm following spontaneous subarachnoid hemorrhage.
- Provides information related to intracranial stenosis or occlusion and cerebral circulatory arrest.
- Gold standard to detect, localize and quantify cerebral embolism in real-time
- Useful for monitoring thrombolysis of acute intracranial occlusions, detecting extracranial ICA stenosis, cerebral microembolism, right-to left cardiac shunts and to assess cerebral vasomotor reactivity
- Mean flow velocities and ratios are the primary criteria for diagnosis of cerebrovascular disease
 - Waveforms obtained at 2-5mm increments with a 3-6mm sample volume size
 - Vessel identification is based on depth and direction of flow
 - Spectral analysis includes peak systolic velocity, end diastolic velocity, systolic upstroke or acceleration time, pulsatility index, and time-averaged maximum mean velocity
 - Calculate mean flow velocity ratios
 - Evaluate Doppler waveform characteristics
 - Document randomly occurring high intensity signals (HITS) associated with embolic phenomenon
- Benefits of TCD when compared to TCDI:

- better beam penetration and vessel detection rates
- can quickly evaluate intracranial arteries to perform provocational testing, monitoring of emboli and vasomotor reactivity testing

REVIEW OF THE DIAGNOSTIC EXAM FINDINGS

The sonographer or examiner should:

- Review data acquired during the exam to ensure a complete and comprehensive evaluation has been performed and documented.
- Explain and document any exceptions or limitations to the protocol.
- Document any changes compared with previous exams.
- Record technical findings on a worksheet or other appropriate method (e.g., computer software) and classify results according to the laboratory diagnostic criteria.
- Document the exam date, clinical indications, performing sonographer and exam summary in patient's medical record.

PRESENTATION OF EXAM FINDINGS

The sonographer or examiner should:

- Provide preliminary results, when necessary, as provided by laboratory specific guidelines.
- Present record of diagnostic images, data, explanations, and technical worksheet to the interpreting physician. Interpretation must be available within two business days.
- Sonographer and interpreting physician name must appear in the final report. Finalized/signed report should be available within four business days.
- Alert vascular laboratory Medical Director or appropriate health care provider when immediate medical attention is indicated based on departmental guidelines and procedures.

EXAM TIME RECOMMENDATIONS

High quality and accurate results are fundamental elements of TCD and TCDI examinations. A combination of indirect and direct exam components is the foundation for maximizing exam quality and accuracy.

- Indirect exam components include:
 - Pre-exam activities: obtaining previous exam data, initiating exam worksheet and paperwork, equipment and exam room preparation, patient assessment and positioning, patient communication
 - Post-exam activities: exam room cleanup, compiling and

processing exam data for preliminary and/or formal interpretation, and exam billing activities.

- Direct exam components include:
 - Equipment optimization and the actual hands-on examination process.
- While study times may vary depending on testing protocols, patient condition, and clinical complexity of the evaluation being performed, these are the times necessary to provide a quality diagnostic evaluation. Listed are the recommended examination times for performing each CPT related to this guideline, which were derived from the direct time inputs from the Resource Based Relative Value Scale (RBRVS)
 - 93886 73 minutes
 - If applicable, +93896 (vasoreactivity study with complete TCD)
 - If applicable, +93897 (emboli detection without IV microbubble, performed with complete TCD)
 - If applicable +93898 (venous-arterial shunt detection with IV microbubble injection)
 - 93888 49 minutes
 - 93890 73 minutes
 - 93892 83 minutes
 - 93893 76 minutes

REFERENCES

Aaaslid, R. (n.d.). *Cerebral hemodynamics and transcranial Doppler training*. <http://www.hemodynamic.com/index.html>

Alexandrov, A. V. (2011). *Cerebrovascular ultrasound in stroke prevention and treatment* (2nd ed.). Wiley-Blackwell.

Alexandrov, A. V., Sloan, M. A., Tegeler, C. H., et al. (2012). Practice standards for transcranial Doppler (TCD) ultrasound: Part II. Clinical indications and expected outcomes. *Journal of Neuroimaging*, 22, 215–224. <https://doi.org/10.1111/j.1552-6569.2010.00523.x>

Alexandrov, A. V., Sloan, M. A., Wong, L. K., et al. (2007). Practice standards for transcranial Doppler ultrasound: Part I—Test performance. *Journal of Neuroimaging*, 17, 11–18. <https://doi.org/10.1111/j.1552-6569.2006.00088.x>

Bulas, D. I., Jones, A. M., Seibert, J. S., et al. (2000). Transcranial Doppler (TCD) screening for stroke prevention in sickle cell anemia: Pitfalls in technique variation. *Pediatric Radiology*, 30(11), 733–738. <https://doi.org/10.1007/s002470000317>

Jawad, N., Kok Hooi, Y., Gulraiz, A., & Jonathan, G. (2013). Transcranial Doppler ultrasound: A review of the physical principles and major applications in critical care. *International Journal of Vascular Medicine*, 2013, Article 629378. <https://doi.org/10.1155/2013/629378>

Jones, A. M., Adams, R. J., Seibert, J. S., Nichols, F. T., et al. (2001). Comparison of transcranial Doppler and transcranial Doppler imaging in patients with sickle cell anemia. *Pediatric Radiology*, 31(7), 461–469. <https://doi.org/10.1007/s002470100427>
Purkayastha, S., & Sorond, F. (2012). Transcranial Doppler ultrasound: Technique and application. *Seminars in Neurology*, 32(4), 411–420. <https://doi.org/10.1055/s-0032-1331812>

Rommer, B., Bellner, J., Kongstad, P., & Sjöholm, H. (1996). Elevated transcranial Doppler flow velocities after severe head injury: Cerebral vasospasm or hyperemia? *Journal of Neurosurgery*, 85(1), 90–97. <https://doi.org/10.3171/jns.1996.85.1.0090>

Sloan, M. A., Alexandrov, A. V., Tegeler, C. H., et al. (2004). Assessment: Transcranial Doppler ultrasonography: Report of the Therapeutics and Technology Assessment Subcommittee of the American Academy of Neurology. *Neurology*, 62, 1468–1481. <https://doi.org/10.1212/WNL.62.9.1468>

Topcuoglu, M. A., Unal, A., & Arsava, E. M. (2010). Advances in transcranial Doppler clinical applications. *Expert Opinion on Medical Diagnostics*, 4(4), 343–358. <https://doi.org/10.1517/17530059.2010.495749>