



VASCULAR TECHNOLOGY PROFESSIONAL PERFORMANCE GUIDELINES

Lower Extremity Arterial Physiologic Evaluations

This Guideline was prepared by the Professional Guidelines Subcommittee of the Society for Vascular Ultrasound (SVU) as a template to aid the vascular technologist/sonographer and other interested parties. It implies a consensus of those substantially concerned with its scope and provisions. The guidelines contain recommendations only and should not be used as a sole basis to make medical practice decisions. This SVU Guideline may be revised or withdrawn at any time. The procedures of SVU require that action be taken to reaffirm, revise, or withdraw this Guideline no later than three years from the date of publication. Suggestions for improvement of this Guideline are welcome and should be sent to the Executive Director of the Society for Vascular Ultrasound. No part of this Guideline may be reproduced in any form, in an electronic retrieval system or otherwise, without the prior written permission of the publisher.

Sponsored and published by:
Society for Vascular Ultrasound
1601 Utica Ave S., Suite 213 Minneapolis, MN
55416
Tel.: 651-288-3431
Fax: 651-290-2266
E-mail: svuinfo@svunet.org
Internet: www.svunet.org

Copyright © by the Society for Vascular Ultrasound, 2026.
ALL RIGHTS RESERVED. PRINTED IN THE UNITED STATES OF AMERICA.
VASCULAR PROFESSIONAL PERFORMANCE GUIDELINE

Updated April 2026

PURPOSE

Segmental pressures, pulse volume recordings, Doppler and photoplethysmography waveform analysis are performed to evaluate lower extremity arteries to determine the presence, severity and location of arterial occlusive disease.

COMMON INDICATIONS

Common indications for performance of arterial segmental pressures and waveforms include, but are not limited to:

- Bruit
- Arterial insufficiency
- Claudication
- Limb / digital ischemia
- Rest pain
- Ulceration
- Peripheral arterial disease (PAD)
- Vasospastic Disease
- Raynaud's Syndrome
- Suspected arterial embolization
- Pre-procedure assessment (e.g., angioplasty, stent, bypass graft, amputation)
- Follow-up assessment of lower extremity revascularization procedures (e.g., angioplasty, stenting, bypass or placement of hemodialysis access conduits)

CONTRAINDICATIONS AND LIMITATIONS

Contraindications for lower extremity arterial physiologic testing include, but are not limited to:

- Morbidly obese patients in whom a high-thigh pressure may not be obtainable
- Casts, dressings, staples, open wounds or ulceration
- Suspected or known acute deep venous thrombosis (DVT)
- Non-compressible vessels with arterial calcification; TBI with digital PPG waveforms must be obtained; the ABI is unreliable in this setting and cannot exclude PAD
- Post-interventional procedures (e.g., stenting, bypass graft)
- Recent surgery or trauma
- Involuntary movement/tremors that may render waveform collection suboptimal or unreliable

PATIENT COMMUNICATION

Prior to beginning the exam, the sonographer or examiner should:

- Introduce self and explain why the examination is being performed and indicate how much time the examination will take.
- Verify the patient's name and date of birth or utilize facility specific patient identifiers.
- Explain the procedure, taking into consideration the age and mental status of the patient.
- Respond to questions and concerns about any aspect of the examination.
- Educate the patient about risk factors and symptoms of peripheral arterial disease.
- Refer specific diagnostic, treatment or prognosis questions to the patient's physician.
- The patient should have rested for at least 15 minutes before beginning the examination.

PATIENT ASSESSMENT

A patient assessment should be completed before the examination. This includes assessment of the patient's ability to tolerate the procedure, an evaluation of any contraindications to the procedure and a physical assessment of bilateral lower extremities, feet and digits. The sonographer or examiner should obtain a complete, pertinent history and review the patient's medical record, if available. A pertinent history includes:

- Previous vascular/cardiovascular surgeries
- Current medications or therapies
- Risk factors for arterial disease
 - Diabetes
 - Hypertension
 - Hyperlipidemia
 - Age
 - Smoking
 - Obesity
 - Cerebrovascular disease
 - Coronary artery disease
 - Family history of PAD
 - Connective tissue disorders
- Presence of symptoms for peripheral arterial disease
 - Claudication
 - Rest pain
 - Ulceration
 - Gangrene
 - Ischemia
 - Hair loss

- Pallor
- Paralysis
- Paresthesia
- Poikilothermia
- Dependent rubor
- Results of other relevant diagnostic procedures
- Review of prior ultrasound examinations

A complete assessment of the patient's history and symptoms should guide the sonographer to:

- Verify that the requested procedure correlates with the patient's clinical presentation
- Determine if supplemental procedures are necessary: exercise treadmill testing, provocative maneuvers or cold challenge
- Perform adjunctive procedures according to the laboratory-specific protocol: auscultation of bruits or palpation of pulses.

PATIENT POSITIONING/PREPARATION

The optimal patient positioning and preparation for performing a lower extremity arterial segmental physiologic evaluation includes the following:

- The extremities positioned at approximately the same level as the heart. The exam is typically performed with the patient in supine position. Position can be altered when performing supplemental procedures
- The examiner should be as close to the extremity as possible to allow for proper ergonomics
- Performed in a warm room to reduce vasoconstriction
- Patient rested for approximately 5-10 minutes prior to beginning examination

INSTRUMENTATION

Non-invasive physiologic studies require separate and distinct equipment from the duplex ultrasound scanner. Instrumentation must allow the display and permanent recording of pressures with bi-directional Doppler analysis of blood flow, plethysmography and/or oxygen tension measurements. The laboratory must be equipped and prepared to routinely acquire and permanently record digital PPG waveforms and toe systolic pressures on every patient.

- Continuous wave (CW) Doppler must provide:
 - Doppler frequencies appropriate for vessels examined
 - Typically range from 4 to 10MHz
 - Direction sensitive Doppler blood flow meter
 - Doppler waveform display with bi-directional flow capabilities
 - Audible output and permanent recording of Doppler waveforms
- Segmental limb plethysmography must provide:
 - Equipment capable of measuring small segmental volume changes and

- making permanent recordings of blood pressure measurements
- Blood pressure cuffs of various sizes for each limb segment and digit evaluated
 - Recommended size is 20% wider than limb diameter
- Pulse volume plethysmography (PVR) must provide:
 - Equipment capable of measuring small limb volume changes and saving permanent PVR waveforms
 - Capability to be calibrated before each exam
 - Blood pressure cuffs of various sizes for each limb segment and digit evaluated
- Photoplethysmography (PPG) must provide:
 - Electrical sensor for signal display
 - Capable of providing a permanent recording of PPG waveforms
- Treadmill exercise/stress testing (if used) must provide:
 - Motorized treadmill capable of maintaining constant speed and incline

EXAMINATION PROTOCOL

Sonographers should follow a standard exam protocol for each lower extremity examination. A complete examination includes bilateral non-invasive physiologic studies (segmental pressures, bidirectional Doppler waveforms, and/or segmental volume plethysmography) of 3 or more levels. A limited bilateral non-invasive physiologic study (segmental pressures, bi-directional Doppler waveforms, and/or segmental volume plethysmography) consist of 1-2 levels.

Throughout each examination, the sonographer should:

- Assess and monitor the patient's physical and mental status, allowing modifications to the procedure plan according to the patient's clinical status.
- Analyze segmental pressure and waveform findings to ensure that sufficient data is provided to the physician to direct patient management and render a final diagnosis.
- Analyze supplemental sonographic findings to ensure that sufficient data is provided to the physician to direct patient management and render a final diagnosis.

Single-Level Lower Extremity Arterial Physiologic Exam

A limited exam includes the measurement of bilateral blood pressures to obtain the ankle brachial index, in combination with either Doppler or plethysmographic waveform analysis at the ankle.

*The **ankle brachial index (ABI)** includes:*

- Measurement of bilateral brachial artery pressures and the higher of the two pressures is used to calculate the ABI

- Measurement of bilateral ankle pressures from the distal posterior tibial (PT) artery and distal anterior tibial (AT)/dorsalis pedis (DP) artery and the higher of the two pressures on each side is used to calculate the ABI
 - The **ankle brachial index** (ABI) is calculated by dividing the highest ankle pressure from each limb by the highest brachial pressure
 - The ABI may be affected by arterial calcification
- Additional information regarding the presence of disease should be obtained by recording toe PPG waveforms and toe pressures
 - The **toe brachial index** (TBI) is calculated by dividing the great toe pressure by the highest brachial pressure

Doppler waveforms are obtained and documented from the:

- Posterior and anterior tibial arteries at the level of the ankle

Plethysmographic waveforms are obtained and documented from the:

- Ankle level
 - PVR waveforms are unaffected by intimal calcification
- Toe pressures, TBI, and digital PPG waveforms are required in patients with diabetes mellitus, chronic kidney disease, known or suspected medial arterial calcification, resting ABI >1.30, non-compressible ankle vessels, rest pain, non-healing wound, ulceration, gangrene, or suspected CLTI. In these patients, an ABI alone, including a normal ABI, does not exclude significant peripheral arterial disease.

Multi-Level Lower Extremity Arterial Physiologic Exam

A complete exam includes the measurement of bilateral blood pressures to obtain the ankle brachial index, in combination with either Doppler or plethysmographic waveform analysis from at least three levels.

The ankle brachial index (ABI) includes:

- Measurement of bilateral brachial artery pressures and the higher of the two pressures used to calculate the ABI
- Measurement of bilateral ankle pressures from the distal posterior tibial (PT) artery and distal anterior tibial (AT)/dorsalis pedis (DP) artery and the higher of the two pressures on each side is used to calculate the ABI
 - The **ankle brachial index** (ABI) is calculated by dividing the highest ankle pressure from each limb by the highest brachial pressure
 - The ABI may be affected by arterial calcification
- Additional information regarding the presence of disease is obtained by recording toe PPG waveforms and toe pressures

- The **toe brachial index** (TBI) is calculated by dividing the great toe pressure by the highest brachial pressure

Segmental pressures at the thigh, calf and ankle levels are obtained in cases with an abnormal ABI:

- Four-cuff technique: includes high thigh, low thigh, calf and ankle pressures
- Three-cuff technique: includes thigh, calf and ankle pressures

Doppler waveforms are obtained and documented from the:

- Common femoral artery
- Mid superficial femoral artery (when indicated)
- Popliteal artery
- Posterior and anterior tibial arteries at the level of the ankle

**Gain settings should be maximized at each level to display Doppler waveform characteristics.*

Plethysmographic waveforms are obtained and documented from the:

- Thigh (includes high thigh and low thigh waveforms if using the four-cuff method)
- Calf
- Ankle
- Toe pressures, TBI, and digital PPG waveforms are required in patients with diabetes mellitus, chronic kidney disease, known or suspected medial arterial calcification, resting ABI >1.30, non-compressible ankle vessels, rest pain, non-healing wound, ulceration, gangrene, or suspected CLTI. In these patients, an ABI alone, including a normal ABI, does not exclude significant peripheral arterial disease.

**Gain settings are optimized to display PVR waveform characteristics, but left unchanged between levels.*

Lower Extremity Arterial Physiologic Exam with Treadmill Exercise/Stress Testing

If indicated and ordered by a qualified physician, a non-invasive physiologic study at rest and following motorized treadmill stress testing may be performed. Treadmill testing is generally set at a constant speed and grade (e.g., 2 mph, 10% grade) for five minutes, or until symptoms occur and the patient is forced to stop.

Resting exam includes the measurement of bilateral blood pressures to obtain the ankle brachial index, in combination with either Doppler or plethysmographic waveform analysis at the ankle. This is followed by treadmill exercise for five minutes, or until symptoms occur and the patient is forced to stop.

Following treadmill exercise, bilateral ankle pressure measurements are obtained at timed intervals with the patient in a supine position:

- Begin the ankle pressure measurements in the symptomatic limb

- or limb with the lowest pre-exercise ankle pressure.
- Use the pedal artery that yielded the higher pre-exercise ankle pressure to record the post-exercise ankle pressure.
- Use the brachial artery that yielded the higher pre-exercise pressure to obtain a post-exercise pressure and calculate the ABI.
- Repeat the post-exercise ABI/pressure measurements at one-to-two-minute intervals for up to 10 minutes, or until ankle pressures return to pre-exercise levels.
 - testing results are based on initial decrease in ankle pressure and the time required for the ankle pressures to return to pre-exercise levels.

Contraindications to treadmill exercise/stress testing

When exercising any patient, the examiner should be familiar with risk factors and contraindications related to this test and aware that the protocol may need to be modified for individual patients. Contraindications for treadmill exercise/stress testing include:

- ABI less than 0.5
- Non-compressible vessels (ABI > 1.3)
- Chest pain (unless physician present)
- Questionable or unstable cardiac status
- Shortness of breath
- Unstable gait
- Hypertension (>180mmHg brachial pressure)

Treadmill exercise testing should be discontinued if:

- Patient terminates the examination due to inability to continue ambulating
- Patient experiences any of the following symptoms:
 - Chest pain; shoulder, neck, jaw, or arm pain
 - Shortness of breath, fatigue, or faintness

REVIEW OF THE DIAGNOSTIC EXAMINATION FINDINGS

The sonographer or examiner should:

- Review data acquired during the lower extremity arterial physiologic exam to ensure that a complete and comprehensive evaluation has been performed and documented.
- Explain and document any exceptions to the protocol (i.e., study omissions or revisions).
- Determine any change in follow-up studies, review previous exam documentation, document any change in status; and, to duplicate prior imaging and Doppler parameters.
- Record technical findings required to complete the final diagnosis on a worksheet or other appropriate method (e.g., computer software), so that

the findings can be classified according to the laboratory diagnostic criteria.

- Document the exam date, clinical indications, sonographer performing the evaluation, and exam summary in the patient's medical record.

PRESENTATION OF EXAMINATION FINDINGS

The sonographer or examiner should:

- Provide preliminary results, when necessary, as provided for by laboratory specific guidelines.
- Present record of diagnostic images, data, explanations, and technical worksheet to the interpreting physician. Interpretation must be available within two business days.
- The sonographer's and interpreting physician's name must appear on the final report. The finalized/signed report should be available within four business days.
- Alert vascular laboratory Medical Director or appropriate healthcare provider when immediate medical attention is indicated based on departmental guidelines and procedures.

EXAMINATION TIME RECOMMENDATIONS

High-quality accurate results are fundamental elements of lower extremity physiologic examinations. A combination of indirect and direct exam components are the foundation for maximizing exam quality and accuracy.

- Indirect exam components include:
 - Pre-exam activities: obtaining previous exam data, initiating exam worksheet and paperwork, equipment and exam room preparation, patient assessment and positioning, patient communication
 - Post-exam activities: exam room cleanup, compiling and processing exam data for preliminary and/or formal interpretation, and exam billing activities
- Direct exam components include:
 - Equipment optimization and the actual hands-on, examination process
- While study times may vary depending on testing protocols, patient condition, and clinical complexity of the evaluation being performed, these are the times necessary to provide a quality diagnostic evaluation. Listed are the recommended examination times for performing each CPT related to this guideline, which were derived from the direct time inputs from the Resource Based Relative Value Scale (RBRVS).
 - 93922 39 minutes
 - 93923 59 minutes
 - 93924 73 minutes

REFERENCES

- Gerhard-Herman, M., [et al.]. (2006). Guidelines for noninvasive vascular laboratory testing: A report from the American Society of Echocardiography and the Society of Vascular Medicine and Biology. *Journal of the American Society of Echocardiography*, 11(8), 955–972. <https://doi.org/10.1177/1358863X06070516>
- Intersocietal Accreditation Commission. (2018). *IAC standards and guidelines for vascular testing*. <https://www.intersocietal.org/vascular/standards/IACVascularTestingStandards2018.pdf>
- Kohler, T. R., & Sumner, D. S. (2010). Vascular laboratory: Arterial physiologic assessment. In J. L. Cronenwett & K. W. Johnston (Eds.), *Rutherford's vascular surgery* (7th ed., pp. 217–233). Saunders Elsevier.
- Kupinski, A. M. (2002). Segmental pressure measurement and plethysmography. *Journal of Vascular Technology*, 26(1), 32–38.
- McPharlin, M. (2012). Lower-extremity arterial air plethysmography evaluation. *Journal of Vascular Ultrasound*, 36(2), 135–142.
- Mohler, E. R., III, Gornik, H. L., Gerhard-Herman, M., Misra, S., Olin, J. W., Zierler, R. E., & [et al.]. (2012). ACCF/ACR/AIUM/ASE/ASN/ICAVL/SCAI/SCCT/SIR/SVM/SVS 2012 appropriate use criteria for peripheral vascular ultrasound and physiological testing part I: Arterial ultrasound and physiological testing: A report of the American College of Cardiology Foundation Appropriate Use Criteria Task Force, American College of Radiology, American Institute of Ultrasound in Medicine, American Society of Echocardiography, American Society of Nephrology, Intersocietal Commission for the Accreditation of Vascular Laboratories, Society for Cardiovascular Angiography and Interventions, Society of Cardiovascular Computed Tomography, Society for Interventional Radiology, Society for Vascular Medicine, and Society for Vascular Surgery. *Journal of Vascular Surgery*, 56(1), e17–e51. <https://doi.org/10.1016/j.jvs.2012.05.054>
- Poe, P. (2012). Lower-extremity arterial continuous-wave Doppler evaluation. *Journal of Vascular Ultrasound*, 36(2), 123–134.
- Sononet. (2017). *Ankle-brachial index: A diagnostic tool for peripheral arterial disease*. https://www.sononet.us/abiscore/documents/ABI_how_to.pdf